



Mats Ljungman, Ph.D.
Department of Radiation Oncology
University of Michigan
ljungman@umich.edu

Award UM1HG009382

Document date: 2020-12-11

Document authors: Michelle Paulsen (tenbroeke@umich.edu), Brian Magnuson (bmagnuso@umich.edu), Ishwarya Venkata Narayanan (ivenkat@umich.edu)

Bru-seq Experiment Protocol

Contents


1. Bromouridine Labeling of Cells
2. Cell Lysis and RNA Extraction
3. Bru-RNA Isolation
4. Stranded Library Preparation

Notes

This protocol is optimized for 150mm dish (1 per sample) with the intent to collect a minimum of 5 million cells (more is better). Bru-RNA makes up ~1% of total RNA, requiring a large amount of starting material for a Bru-RNA library. For BruUV-seq and BruChase-seq, the % Bru-RNA is expected to be even less.

The main differences in protocols between Bru-seq and the derivative techniques BruUV-seq and BruChase-seq are found in the labeling protocol. Bru-seq with IR treatment protocol is identical to Bru-seq (no treatment) beginning with the bromouridine labeling section.

Visual Protocol



IR Treatment

Equipment

IC-320 Biological Irradiator (Kimtron Inc., Oxford, CT)

Steps

1. Following standard safety and operating practices for the irradiator, dose cells with 5 Gy x-rays (66 sec at dose rate of 4.5223 Gy/min).
2. Incubate at 37°C for 1 hour.
3. Immediately label cells with bromouridine (next section).

Bromouridine Labeling of Cells

Reagents

50mM Bromouridine (Aldrich; cat # 850187) in PBS

Steps

1. Add Bru.
 - a. Adherent cells: Remove 10ml of media from each plate of cells to a sterile tube and add BrU to a final concentration of 2mM. Discard remaining media. Add back BrU-containing media to the plate.
 - b. Suspension cells: Measure media volume and add Bru to a final concentration of 2mM.
2. Incubate at 37°C for 30 minutes.
3. Lyse cells immediately (next section).

Cell Lysis and RNA Extraction

Reagents/Materials

Round bottom polypropylene test tube (e.g. Falcon 352059)

TRIzol (Invitrogen; cat # 15596-026)

Chloroform

Isopropanol

75% Ethanol in DEPC-H₂O

DEPC water

Turbo DNA free kit (ThermoFisher; cat #AM1907)

1. Lyse cells in TRIzol.
 - a. Adherent cells: Aspirate media and immediately add 3ml TRIzol to the plate to lyse cells. Scrape cells and transfer lysate to a round bottom tube.
 - b. Suspension cells: Transfer cells and media to a round bottom tube, spin down cells, aspirate media, and resuspend pellet in 3ml TRIzol.
2. Store at -80°C.
3. Thaw to room temperature.
4. Follow TRIzol protocol for RNA extraction.
5. Resuspend pellet in 100µl DEPC-water and incubate at 55°C for 10min to fully dissolve RNA.
6. Follow Turbo DNA free kit protocol.
7. Store total RNA at -80°C.

Bru-RNA Isolation

Reagents/Materials

Goat anti-mouse Dynabeads (Invitrogen cat# 110.33)

Magnetic stand for Eppendorf tubes

Anti-BrdU (BDPharmingen 555627)

0.1% BSA in DEPC-PBS

RNaseOUT Ribonuclease inhibitor (Invitrogen; Cat #10777019) : Diluted to 0.5x using 0.1% BSA in DEPC-PBS

RNA spike-in cocktail (see spike-in reference documentation)

1. Use 50µl aliquot of magnetic beads per sample. Capture beads with magnet and aspirate storage buffer.

2. Wash 3x with 200µl 0.1% BSA in PBS, capturing beads between each wash.
3. Resuspend in 200µl 0.1% BSA in PBS
4. Add 2.5µg (5µl) anti-BrdU monoclonal antibody to beads, then add 1.0µl diluted RNaseOUT.
5. Conjugate by incubating at room temperature with gentle rotation for 1 hour.
6. Wash beads 3x with 200µl 0.1% BSA in PBS, capturing beads between each wash.
7. Resuspend beads in 200µl 0.1% BSA in PBS plus 1.0µl diluted RNaseOUT.
8. Thaw aliquot of RNA spike-in cocktail (on ice) and select appropriate volume to add per 100µg of sample total RNA:
 - a. For Bru-seq, use 25µl per 100µg of sample total RNA.
 - b. For BruUV-seq, use 20% of the Bru-seq amount (5µl).
 - c. For BruChase-seq (6h chase), use 25% of the Bru-seq amount (6.25µl).
 - d. For BruChase-seq (2h chase), use the same amount as for Bru-seq (25ul).
9. Combine total RNA, spike-ins, and water to a total volume of ~200µl.
10. Incubate isolated RNA to 80°C for 10 minutes.
11. Add RNA solution to conjugated beads, incubate at room temp with gentle rotation for 1 hour.
12. Wash beads with 400µl 0.1% BSA in PBS, rotating for 1 minute.
13. Briefly centrifuge, capture beads, and discard supernatant.
14. Add 200µl 0.1% BSA in PBS, flick to mix. Capture beads, discard supernatant.
15. Add 200µl DEPC-PBS, flick to mix. Capture beads, discard supernatant.
16. Add 35µl DEPC-H2O to beads. Use less for a more concentrated sample. Immediately transfer to a new (siliconized or LoBind) tube.
17. Incubate at 96°C for 10 min in a thermomixer set to 550rpm. Briefly centrifuge, then capture beads and keep the supernatant, which contains the Bru-RNA.
18. Measure RNA concentration and store at -80°C.

Stranded Library Preparation

Reagents/Materials

200ng Bru-RNA in 16µl (can use more)
 5x First-Strand buffer
 Actinomycin D
 DTT
 Random Primers (3ug/µl)
 Fast Select (Qiagen; cat# 333180)
 Superscript II (Invitrogen; cat # 18064-014)
 Random primer (Invitrogen; cat # 48190-011)
 NEBuffer 2 (NEB; cat # B7002S)
 RNase H (Invitrogen; cat # 18021-071)
 dUTP (Roche; cat # 11934554001)

DNA polymerase 1 (Invitrogen; cat # 18010-025)
T4 Ligase Buffer (NEB; cat # B0202S)
T4 DNA polymerase (NEB; cat # M0203S)
T4 PNK (NEB; cat # M0201S)
Klenow 3' - 5' exo (NEB; cat # M0212S)
Quick Ligation Kit (NEB; cat # M2200S)
NuSieve 3:1 agarose (Lonza; cat # 50090)
Accugene 10X TAE (Lonza; cat # 50844)
Qiaex II Kit (Qiagen; cat # 20021)
Phusion Master Mix (ThermoFisher; cat # F548S)
USER (NEB; cat # M5505L)
RNA Clean beads (Fisher; cat # NC0068576)
AmPure beads (Fisher; cat # NC9933872)
dNTP set (Invitrogen; cat # 10297-018)
RNaseOUT (Invitrogen; cat # 10777-019)
100bp DNA ladder (NEB; cat # N3231S)

Notes

Steps denoted with *** are stopping points where samples can be stored at -20°C for up to seven days.

1. Remove rRNA with FastSelect.
 - a. Prepare pre-mix: 8µl First-Strand buffer, 1µl Random Primers, 0.10µl Fast Select (per sample).
 - b. Add 9.10 µl of pre-mix, then 16µl RNA, to each PCR tube. Mix well.
 - c. Run FastSelect protocol in thermal cycler. Cool to 4°C.
2. Synthesize First Strand cDNA.
 - a. Mix (per sample): 6.9µl ddH₂O, 4.0µl 100mM DTT, 0.8µl 25mM dNTP, 0.8µl ActinomycinD (2.5µg/µl), 0.5µl RNaseOUT, 2.0µl Superscript II
 - b. Run thermal cycler program: 1) 25°C for 10 minutes, 2) 42°C for 50 minutes, 3) 70°C for 15 minutes, 4) hold at 4°C.
 - c. Clean up with AMPure RNAClean beads.
 - i. Add 72µl AMPure beads in a microcentrifuge tube, then add cDNA mixture.
 - ii. Bind at 25°C for at least 10min (with shaking) and capture beads.
 - iii. Remove and discard supernatant from each sample using micropipettor.
 - iv. Without disturbing the beads, add ~200µl freshly prepared 80% ethanol.
 - v. Incubate at room temperature for 30sec, then discard supernatant.
 - vi. Repeat ethanol wash.
 - vii. Spin briefly and remove any remaining ethanol.
 - viii. Dry beads at 37°C, until just cracked.

- d. Resuspend the dried pellet with 32µl of 5mM Tris pH 8.0. Mix well. Incubate at 28°C for 6min (with shaking).
- e. Capture beads and transfer 30µl of the supernatant to a new PCR tube.
3. Synthesize Second Strand cDNA.
 - a. Mix (per sample): 26.48µl ddH₂O, 15µl NEBuffer 2, 0.9µl 25mM dG+dA+dU+dC mix, 0.75µl RNase H, 1.875µl DNA polymerase
 - b. Add 45µl mixture to each 30µl First Strand product (from Step 2) and mix well.
 - c. Incubate using thermal cycler at 16°C for 2 hours.
 - d. Clean up with 112.5µl AMPure beads. Proceed as per step 2.c. above.
 - e. Resuspend the dried pellet with 27µl 5mM TRis pH 8.0. Mix well. Incubate at 28°C for 6min (with shaking).
 - f. Capture beads and transfer 25µl of the supernatant (now ds cDNA) to a new microcentrifuge tube.***
4. Perform End Repair.
 - a. Mix (per sample): 14.7µl ddH₂O, 5µl T4 DNA Ligase Buffer, 0.8µl 25mM dNTPs, 2µl T4 DNA Polymerase, 0.5µl Klenow, 2µl T4 PNK
 - b. Add 25µl mixture to each 25µl Second Strand product (from Step 3) and mix well.
 - c. Incubate samples at 20°C for 30 minutes.
 - d. Clean up with 75µl AMPure beads. Proceed as per step 2.c. above.
 - e. Resuspend the dried pellet in 32 µl 5mM Tris pH 8.0. Incubate samples at 28°C for 6min (with shaking).
 - f. Capture beads and transfer 30µl of the supernatant to a new microcentrifuge tube.***
5. Adenylate 3' Ends.
 - a. Mix (per sample): 11.5µl ddH₂O, 5µl NEBuffer 2, 1µl 10mM dATP, 2.5µl Klenow Exo(-)polymerase
 - b. Add 20µl mixture to 30µl eluted cDNA (product of Step 4) and mix well.
 - c. Incubate samples in a thermomixer at 37°C for 30min.
 - d. Clean up with 75µl AMPure beads. Proceed as per step 2.c. above.
 - e. Resuspend the dried pellet in 20 µl 5mM Tris pH 8.0. Incubate samples at 28°C for 6min (with shaking).
 - f. Capture beads and transfer 18µl of the supernatant to a new microcentrifuge tube.***
6. Ligate Adapters.
 - a. Anneal ligation adapters: thaw adapter aliquot, incubate at 55°C for 1min in a thermomixer, and cool to 22°C.
 - b. Mix (per sample): 20µl 2x Quick Ligase Buffer, 1µl Quick ligase
 - c. Add 1.5µl annealed ligation adapter to 18µl eluted cDNA (product of step 5), then add 20µl ligase mix. Mix well.
 - d. Incubate at room temp for 30 minutes.
 - e. Clean up with 40µl AMPure beads. Proceed as per step 2.c. above.
 - f. Resuspend the dried pellet in 32 µl 5mM Tris pH 8.0. Incubate samples at 28°C for 6min (with shaking).

- g. Capture beads and transfer 30µl of the supernatant to a new microcentrifuge tube.
 - h. Add 5µl gel-loading buffer.***
 7. Perform size selection by Agarose gel (3% gel using NuSieve 3:1 agarose).
 - a. Remove buffer from wells before loading.
 - b. Run gel in TAE (do not cover gel with buffer) at 75V for 80 min.
 - c. Excise gel slices at 500bp. Reserve a larger slice as a backup.
 - d. Purify gel slices using the Qiagen QIAEXII kit. Elute in 42µl resuspension buffer to recover 40µl. Reserve 20µl as a backup.
 8. Perform uridine digestion and enrich DNA fragments.
 - a. Add 1µl USER to each 20µl sample. Incubate in PCR machine 37°C for 15min, then cool to 4°C.
 - b. Prepare MasterMix/TIPP: Add 1µl TIPP per 250µl MasterMix.
 - c. Put samples on ice.
 - d. Add 1.5µl adapter primer mix.
 - e. Add 25µl Phusion Master Mix/TIPP and mix well.
 - f. Incubate samples (hot start) in a thermal cycler using the program:
 - i. 98°C for 30sec
 - ii. 11-14 cycles of:
 1. 98°C for 10sec
 2. 60°C for 30sec
 3. 72°C for 30sec
 - iii. 72°C for 5min
 - iv. Hold at 10°C
 - g. Clean up with 38µl AMPure beads. Proceed as per step 2.c. above.
 - h. Resuspend the dried pellet in 36µl 5mM Tris pH 8.0. Incubate samples at 28°C for 6min (with shaking).
 - i. Capture beads and transfer 36µl of the supernatant to a new low-bind microcentrifuge tube.
 9. Measure DNA concentration (Nanodrop or similar).
 10. Run 3µl of each library on a confirmation gel (1.5% agarose).