

IGSB transfection of K562 cell-Lipofectamine 3000

1. Seed 500,000 pure K562 cells in 12 well plate, in 1mL complete RPMI media. Include a control.
2. Label a 1.5mL centrifuge tube for each transfection and add 50 μ L opti-MEM media to each.
3. Add 1 μ g of BAC DNA to each tube and mix by flicking tube
4. Mix P3000 reagent gently, add 2 μ L to each 1.5mL tube. Incubate at RT for 5 min.
5. Add 3 μ L Lipofectamine 3000 to a new 1.5 mL centrifuge tube containing 50 μ L opti-MEM. Mix by flicking.
6. Add diluted DNA to diluted Lipofectamine 3000 and incubate for 5 min at RT.
7. Add entire volume of transfection mixture drop-wise to appropriate cells in the 12 well plate.
8. **Day after transfection:** Spin down the cells at 500 x g for 5 min. at room temperature. Resuspend the cells in 2mL complete media and place back into the 12 well plate.
9. **Two days after transfection:** Spin down the cells again at 500 x g for 5 min at room temp. and place in 6 well plates with 4mL of complete medium per well. Add Geneticin G418 to a concentration of 50 μ g/mL to begin light selection.
10. **Five days after transfection:** Replace the media with full selection media (500 μ g/mL G418) and place back into 6 well plate.
11. Continue to replace selection media every two days (usually M-W-F).
12. **Nineteen days after transfection:** expand the cells to a T-25 flask and bring volume up to 10mL.
13. When negative control has completely died (~2.3 to 3 weeks), expand cells to a T-75 flask.
14. Once cells have reached log growth phase, doubling approximately every day. Collect 7x10⁶ cells and centrifuge 500 x g for 5 min.
15. Resuspend the pellets in 7mL recovery cell culture medium with 5% DMSO and aliquot into 7 cryogenic vials.
16. Place tubes in "Mr. Frosty" Nalgene box with isopropanol and place at -80°C. Once frozen, the cells can be placed in liquid nitrogen for storage.
17. **Crosslinked Cells:** With the cells suspended in RPMI media, add formaldehyde to a final concentration of 1%. Incubate cells for 10 min.
18. Add a 1:20 volume of 2.5M Glycine, and keep on ice.
19. Wash twice with cold PBS, centrifuged as above. Pellets can be snap frozen and stored at -80°C.

IF protocol

1. Place 200 μ L of cells in a 1.5mL tube.
2. Centrifuge at 2000 rpm and wash twice with 500 μ L of PBS with 4% formaldehyde.
3. Incubate for 20 min at RT.
4. Wash twice with 500 μ L of PBS with 0.1% BSA.
5. Centrifuge and replace wash buffer with 500 μ L Blocking Buffer.
6. Incubate for 45 min at RT.
7. Remove blocking buffer and add 500 μ L Dilution Buffer with a 1:1000 dilution of Ab290 (Abcam). Stock can be stored at 4°C for several months or 5 times.
8. Incubate overnight at 4°C with rotation.
9. Centrifuge and wash twice with 500 μ L Wash buffer.
10. Resuspend in 500 μ L Dilution buffer with 1 μ L of Alexa488. Incubate for 1 hr, rotating in the dark.
11. Centrifuge and wash twice with 500 μ L Wash buffer.
12. Leaving ~30 μ L of Wash buffer behind add two drops Prolong mounting media with DAPI (Life Tech. P36935) and transfer to a slide.
13. Place coverslip and allow to dry at RT in the dark for several hours. Store at 4°C.

Blocking Buffer:

PBS with 0.3 % Triton X-100, 10% Goat serum.

Wash Buffer:

PBS with 0.1% BSA.

Dilution Buffer:

PBS with 1% BSA, 1% goat serum, 0.3% Triton X-100.