

SOP: Propagation of NCI-H226 Human Lung Mesothelioma Cells
Date modified: 7/13/2012
Modified by: T. Canfield (UW)

Ordering Information

NCI-H226 Human Lung Mesothelioma Cells can be ordered from ATCC as a frozen ampoule with 0.9×10^6 cells per 1mL volume. This is an adherent metastatic squamous cell carcinoma epithelial cell line from a human male.

Name: NCI-H226—Human Lung Mesothelioma Cells
ATCC #: CRL-5826

Materials List

1. RPMI-1640 Medium (ATCC, Cat# 30-2001)
2. Characterized Fetal Bovine Serum (HyClone, Cat# SH30071)
3. Penicillin-Streptomycin Solution (200X) (Cellgro, Cat# 30-001-CI)
4. T25, T75, T225 tissue culture flasks
5. Corning conical centrifuge tubes (15mL and 50mL)
6. Graduated pipets (1, 5, 10, 25, 50mL)
7. Phosphate Buffered Saline (1X PBS) (Cellgro, Cat# 21-040-CM)
8. Accutase Enzyme Cell Detachment Medium (EBiosciences, Cat# 00-4555)
9. Freezing Medium (Growth medium containing 5% DMSO)
10. DMSO, Hybri-Max (Sigma-Aldrich, Cat# D2650)
11. Cryovials (Nunc, Cat# 368632)
12. Cryo 1°C Freezing Container (Nalgene Cat# 5100-0001)
13. Eppendorf Centrifuge 5810R
14. Revco UltimaII -80°C Freezer
15. Thermolyne Locator 4 Liquid Nitrogen Freezer
16. Hemocytometer
17. Micropipet w/ P20 tips
18. Microscope

Growth Medium for NCI-H226

RPMI-1640 Medium
10% Characterized FBS
Pen-Strep (1X)

Procedure

A. Receipt of Frozen Cells and Starting Cell Culture

- 1) Immediately place frozen cells in liquid nitrogen storage until ready to culture.
- 2) When ready to start cell culture, quickly thaw ampoule in a 37°C water bath.
- 3) As soon as ice crystals disappear, swab outside surface of the ampoule with 70% ethanol, then dispense contents of ampoule into a 15mL Corning centrifuge tube containing 9mL complete culture medium.

- 4) Pellet cells at 125 x g for 7 minutes (4°C).
- 5) Resuspend cell pellet in 10mL complete culture medium and dispense into a T25 flask.
- 6) To incubate the culture, place the flask in a 37°C, 5% CO₂ humidified incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Aspirate medium.
- 3) Wash cells with room temperature 1X PBS.
- 4) Add 5mL of Accutase and return to incubator for 5-10 minutes, or until cells detach.
- 5) Immediately remove cells, rinse flask with room temperature 1X PBS to collect residual cells, and pellet at 300 x g for 5 minutes (4°C).
- 6) Gently re-suspend cell pellet in prewarmed medium.
- 7) Perform 1:4 to 1:8 cell split as needed.
- 8) Record each subculture event as a passage.
- 9) Note: use increasing amounts of Accutase with increasing sizes of tissue culture vessel to sub-culture.

C. Maintenance and Generation of Seed Stocks

- 1) Change media the day after seeding and every 2-3 days thereafter. Use 10mL of growth medium per T25 flask, 20mL of growth medium per T75 flask, and 50mL of growth medium per T225 flask.
- 2) Following first or second passage after receipt of cells and with sufficient number of cells to continue maintenance and expansion, the major portion of the flasks should be sub-cultured using Accutase as above under “Sub-culture” and a small portion should be set aside as a seed stock. The cell pellet for the seed stock should be resuspended in freezing medium.
- 3) Cells in freezing medium are dispensed into cryovials (2 million cells per 1 mL aliquot) and frozen at -80°C in a Nalgene Cryo 1°C freezing container overnight.
- 4) Cryovials are transferred the next day to liquid nitrogen freezer for long-term storage.

D. Harvest

- 1) Passage cells until the desired number of cells is reached.
- 2) Remove cells from flasks according to protocol described above under “Sub-culture”, using 15mL of Accutase per T225 flask.
- 3) Examine viability using Trypan blue staining (SOP TP-7).