



# Transfection and induction of FLAG-HA tagged construct into cell culture:SC:Graveley

Revision as of 18:07, 13 June 2012 by [N.Washington](#) (Talk | contribs | block) (diff) ← Older revision | Current revision (diff) | Newer revision → (diff)

## navigation

- [Main Page](#)
- [modENCODE Home](#)
- [Recent changes](#)
- [Help](#)

## dcc

- [Design Documents](#)
- [Meetings](#)
- [Resources](#)
- [Progress Reports](#)
- [Developer Tips](#)

## modencode community

- [Consortium Meeting 2011](#)
- [Projects](#)
- [Meetings](#)
- [Publications](#)
- [Working Groups](#)
- [Add a user](#)
- [References](#)
- [Policies/Standards](#)
- [Submission instructions](#)
- [Md Course Review](#)
- [Photos](#)

## analysis wgs

- [modENCODEAWG](#)
- [ENCODEAWG](#)
- [pw="encode/human"](#)

## protocols

- [Main List](#)
- [Submit/Edit](#)
- [Help](#)

## reagents

- [Main List](#)
- [Submit/Edit](#)
- [Help](#)

## experiment descriptions

- [Main List](#)
- [Submit/Edit](#)
- [Help](#)

## search




## toolbox

- [What links here](#)
- [Related changes](#)
- [Upload file](#)
- [Special pages](#)
- [Printable version](#)
- [Permanent link](#)

## Protocol Text

[\[edit\]](#) [\[rich edit\]](#)

FLAG-HA constructs were prepared according to [Construction of FLAG-HA expression clones:SC:Graveley](#) for a specific target. Plasmid DNA was then prepared using the PureLink™ HQ Mini Plasmid Purification Kit [Invitrogen K2100-01], and eluted with 50ul sterile water into a 1.5ml DNA LoBind microfuge tube (Eppendorf). The DNAs were quantified and assayed for A260/A280 ratio >1.8 using a NanoDrop spectrophotometer (ThermoFisher). Drosophila S2R+ cells were grown at 25°C in non-tissue culture treated polystyrene flasks (Corning Incorporated). Each clone was transiently transfected into a single 54ml culture of Drosophila S2R+ cells (1 x 10<sup>6</sup> cells per ml) grown in a T-150 flask in Schneider's Drosophila Media (GIBCO 21720) with 10% heat-inactivated fetal bovine serum (GIBCO 26140079). Twelve micrograms of each DNA was combined with 300ul of Effectene (Qiagen) following the manufacturer's protocol. Twenty four hours after transfection, expression of the tagged protein was induced with 0.35mM CuSO<sub>4</sub>. This level of CuSO<sub>4</sub> has been tested to induce a low-to-medium level of recombinant protein expression for a majority of representative clones. Cells were grown for 24 hours after induction.

## Validation Form

[\[edit\]](#) [\[rich edit\]](#)

(This section to be completed by Project Bioinformatics contact. Toggle the 'help' link below, or contact your DCC Liaison with questions.)

## Protocol "Transfection and induction of FLAG-HA tagged construct into cell culture:SC:Graveley" (Version 4)

Protocol Type:  ?

Input type:  ?

Ouput type:  ? [obi-biomaterial:cell cult...](#)

Short Description:  ?

Each clone was transiently transfected into a single 54ml culture of Drosophila S2R+ cells (1 x 10<sup>6</sup> cells per ml) grown in a T-150 flask in Schneider's Drosophila Media (GIBCO 21720) with 10% heat-inactivated fetal bovine serum (GIBCO 26140079). Twelve micrograms of each DNA was combined with 300ul of Effectene (Qiagen) following the manufacturer's protocol. Twenty four hours after transfection, expression of the tagged protein was induced with 0.35mM CuSO<sub>4</sub>. This level of CuSO<sub>4</sub> has been tested to induce a low-to-medium level of recombinant protein expression for a majority of representative clones. Cells were grown for 24 hours after induction.

URL:  ?

References:  ?

Species:  ?

Lab:  ?

[Autogenerated minimum SDRF template](#)

Please use this page's permanent link when referencing it in data submission (e.g. in the IDF):

[http://wiki.modencode.org/project/index.php?title=Transfection\\_and\\_induction\\_of\\_FLAG-HA\\_tagged\\_construct\\_into\\_cell\\_culture:SC:Graveley&oldid=78429](http://wiki.modencode.org/project/index.php?title=Transfection_and_induction_of_FLAG-HA_tagged_construct_into_cell_culture:SC:Graveley&oldid=78429)

IE Users: Right-click and choose 'Copy Shortcut' to copy the permalnk URL to the clipboard.

Categories: [Protocol](#) | [Celniker Group Protocol](#)