

HARVEST and RNA ISOLATION

Name _____

Date _____

Check One:

- wells 1-6 (hrs 1-6) wells 1-6 (30min, 7, 8, 10, 12hrs)

Consumables:

RNeasy mini kit (Qiagen, cat# 74104)

STEP 1

- Remove RNeasy mini kit (Qiagen, cat# 74104)
- Ensure station is clean: pipettes, bench, and gloves sprayed with RNase free then EtOH, tubes are labelled
- Add 10uL β -mercaptoethanol to 1mL Buffer RLT. The mixture can be stored @ RT for 1 month
**For 6 samples add 30uL β -ME to 3mL Buffer RLT (only 2.1mL needed)
- Wash cells twice with 3mL PBS
- add 350uL Buffer RLT with β -ME into wells
- use a syringe to pipette the cells up and down at least 5 times while rinsing the wells
- transfer cells by pipetting to rnase/dnase free tubes on ice
**SAFE STOPPING POINT - snap freeze cells and stored @ -80C

Name _____

Date _____

Check One:

- wells 1-6 (hrs 1-6) wells 1-6 (30min, 7, 8, 10, 12hrs)

Consumables:

RNase free DNase set (Qiagen cat#79254)

RNeasy mini kit (Qiagen, cat# 74104)

STEP 2

- **If cells were frozen, thaw on ice before proceeding to next step
- **remove RNase free DNase set, keep on ice.
- **If not already done, add 4 volumes of ethanol (96–100%) to RPE as indicated on the bottle to obtain a working solution
- **If not already done dissolve the lyophilized DNase I (1500 Kunitz) in 550 µl of the RNase- free water. (If new, remove the silver tab and use a syringe to add water). Mix gently by inverting the vial. Do not vortex. Keep on ICE.

- Mix RDD and DNaseI as follows (Do NOT vortex):
*Freeze the RNase free DNase set after reconstituting the DNaseI

	vol/sample (uL)	6	# of samples (5% excess built in)
RDD	70	441	
DNaseI	10	63	
total	80	510	

- Add 350uL of 70% EtOH to each sample. Pipette up and down to mix
- Transfer 700uL of sample to a Rneasy mini spin column placed in a 2mL collection tube
- Centrifuge for 15s @10,000 RPM. Discard flow-through
- Add 350uL RW1 to the RNeasy spin column. Centrifuge for 15s @10,000 RPM. Discard flow-through
- Pipette 80uL RDD/DNaseI mix to each sample, be sure to cover but not disturb the membrane
- Incubate on bench for 15min@RT
- After 15min, add 350uL RW1. Spin 15sec @10,000 RPM. Discard flow-through
- Add 500uL buffer RPE to spin column. Centrifuge for 15s @10,000 RPM. Discard flow-through
- Add 500uL buffer RPE to spin column. Centrifuge for 1min @10,000 RPM. Discard flow-through
- Centrifuge at full speed for and additional 2min to dry the membrane
- Place the RNeasy spin column in a new 1.5mL collection tube
- Add 30uL Rnase free water to the membrane without disturbing the membrane
- Allow tubes to sit with lids closed @RT for 1min.
- After 1min. Spin tubes for 1min @ 15,000RPM
- Place tubes on ice. Check on Nanodrop.
- *if not proceeding to next step, samples can be stored at -80C
- Create sample plate map

- Plate out 20ng/uL in a volume of 55uL on ICE (only 1ug in 50uL is needed for lib prep)
- Check quality of RNA using Agilent RNA screen tape