

SOP: Propagation of ELF-1 Human Embryonic Stem Cell Line
Date modified: 04/02/2013
Modified by: J. Hesson (Carol Ware Lab, UW)

Cell Line Information

Human embryonic stem cell line ELF-1 was derived from a human embryo in the Ellison Stem Cell Core at the University of Washington, Seattle, WA. [NIH Stem Cell Registry (0156); will be banked at WiCell for distribution].

Materials List

1. DMEM/F12 Media (1X), liquid 1:1 (Invitrogen, Cat# 11320-082)
2. Knockout Serum Replacer (Invitrogen, Cat# 10828-028)
3. MEM Sodium Pyruvate Solution, 100X (Invitrogen, Cat# 11360-070)
4. Non-essential Amino Acids, 100X (Invitrogen, Cat# 11140-050)
5. Penicillin/Streptomycin, liquid (Invitrogen, Cat# 15070-063)
6. 2-Mercaptoethanol (Sigma-Aldrich, Cat# M7522)
7. Dulbecco's Modified Eagle Medium (DMEM) (1X), liquid (high glucose), GlutaMAX™ (Invitrogen, Cat# 10566-024)
8. Fetal Bovine Serum (Atlanta Biologicals, Cat# S10250)
9. Gelatin from porcine skin Type A, powder, cell culture tested (Sigma-Aldrich, Cat# G1890)
10. BD Matrigel™ Basement Membrane Matrix, Growth Factor Reduced (GFR) (BD, Cat# 354230)
11. Human LIF (Millipore, Cat# LIF1010)
12. MEK Inhibitor (Selleck Chemicals, Cat# S1036)
13. GSK3 Inhibitor (Selleck Chemicals, Cat# S1263)
14. 35mm and 10cm Tissue Culture Dishes
15. Conical Polypropylene Centrifuge Tubes (15mL and 50mL)
16. Graduated Serological Pipets (1, 5, 10, 25, 50mL)
17. Dulbecco's Phosphate-Buffered Saline (D-PBS) 1X, liquid (Invitrogen, Cat# 14190-250)
18. 0.05% Trypsin-EDTA (1X), Phenol Red (Invitrogen, Cat# 25300-120)
19. Accutase Enzyme Cell Detachment Medium (EBiosciences, Cat# 00-4555)
20. cOmplete™, EDTA-free, Protease Inhibitor Cocktail Tablets (Roche Diagnostics Corp., Cat# 05056489001); use 1 tablet per 50mL solution.
21. Freezing Medium (60% Growth Medium, 30% Knockout Serum Replacer, 10% DMSO)
22. DMSO, Hybri-Max (Sigma-Aldrich, Cat# D2650)
23. Freezing Straws (Veterinary Concepts, Cat# 04170)
24. Gamma-irradiated MEF Cells (generated in the Ellison Stem Cell Core, and frozen)
25. Beckman Coulter Allegra X-22 Centrifuge
26. Bio-Cool Controlled Rate Freezer
27. New Brunswick Cell Counter
28. Scissors
29. Microscope

Growth Medium for ELF-1 Cells

500mL DMEM/F12 (1:1) Medium
100mL Knockout Serum Replacer
6mL Pen/Strep
6mL MEM Sodium Pyruvate Solution
6mL NEAA (Non-essential Amino Acids)
0.6mL 0.1M 2-Mercaptoethanol
0.8 μ M MEK Inhibitor
1.5 μ M GSK3 Inhibitor
10 μ g/mL Human LIF

Freezing Medium for ELF-1 Cells

3mL Growth Medium
1.5mL Knockout Serum Replacer
0.5mL DMSO

Plating Medium for Irradiated MEF Cells

500mL DMEM GlutaMAX™ (High Glucose) Medium
60mL Fetal Bovine Serum
6mL Pen/Strep

Procedure

A. Thawing Frozen Cells and Starting Cell Culture

- 1) When ready to start cell culture, quickly thaw 1 straw in a small container of room temperature water.
- 2) After about 10 seconds, spray straw and a pair of scissors with 70% ethanol, then dispense contents of straw by cutting the top and bottom off and allowing contents to empty into a 15mL conical centrifuge tube containing 8mL complete culture medium.
- 3) Pellet cells at 300 x g for 3 minutes.
- 4) Re-suspend cell pellet in 4mL complete culture medium and dispense into two 35mm tissue culture dishes with 1.5×10^5 irradiated MEFs (9.0×10^5 irradiated MEFs on 10cm tissue culture dish), plated overnight on 0.1% gelatin-coated tissue culture dishes.
- 5) To culture, place the dish in a 37°C, 5% O₂, 5% CO₂, 90% N₂, humidified incubator.

B. Sub-culture

- 1) Propagate cells for 3-4 days, changing medium every 1-2 days.
- 2) Aspirate medium.
- 3) Wash cells with 1X PBS.
- 4) Add 1.5mL of Trypsin to 35mm tissue culture dish (3mL for 10cm tissue culture dish) and let sit at room temperature for 5-10 minutes, or until cells detach.

- 5) Immediately remove cells, rinse tissue culture plate with equal amount of growth medium to collect residual cells, and pellet at 300 x g for 3 minutes.
- 6) Gently re-suspend cell pellet in growth medium.
- 7) Perform 1:5 to 1:10 split every 3-4 days.

C. Harvest

- 1) Passage cells until the desired number of cells is reached. On the last passage before collection, plate cells on Matrigel (diluted 1:30 in PBS) coated tissue culture dishes.
- 2) At time of harvest, rinse plates with PBS.
- 3) Add 3mL Accutase and incubate for 10-15 minutes at 37°C.
- 4) Remove cells to 15mL conical centrifuge tube and rinse dish with culture medium to collect residual cells.
- 5) Pellet cells at 300 x g for 5 minutes.
- 6) Wash pellet in PBS supplemented with cOmplete™, EDTA-free, Protease Inhibitor Cocktail Tablet (1 tablet per 50mL PBS).
- 7) Count number of cells and proceed to **SOP for cultured cells: nuclei, DNaseI treatment, crosslinking, and preserving cells for RNA** and **SOP: Cell purification using Percoll step gradients**.