

SOP: Namalwa cell culture and infection by Sendai Virus
Date modified: 03/25/2013
Modified by: Junqiang Ye (Maniatis Lab, Columbia University)

Ordering Information

Human B cell line Namalwa can be ordered from ATCC. This is a suspension cell line.

Name: Namalwa—Human B Lymphocyte Cells
ATCC #: CRL-1432

Materials List

1. RPMI-1640 medium with 2mM L-glutamine adjusted to contain 1.5g/L sodium bicarbonate, 4.5g/L glucose, 10mM HEPES, and 1.0mM sodium pyruvate, (ATCC, Cat# 30-2001)
2. Characterized Fetal Bovine Serum (HyClone, Cat# SH30071)
3. Penicillin-Streptomycin Solution (200X) (Corning/Cellgro, Cat# 30-001-CI)
4. Sendai Virus Cantell strain (Charles River Laboratories, Cat# 10100772, viral titer > 2000 HAU/mL)
5. T75 tissue culture flasks
6. Corning conical centrifuge tubes (15mL and 50mL)
7. Graduated pipets (1, 5, 10, 25, 50mL)
8. Phosphate Buffered Saline (1X PBS) (Corning/Cellgro, Cat# 21-040-CM)
9. Freezing Medium (Growth medium containing 10% DMSO)
10. DMSO, Hybri-Max (Sigma-Aldrich, Cat# D2650)
11. Cryovials (Nunc, Cat# 368632)
12. Cryo 1°C Freezing Container (Nalgene Cat# 5100-0001)
13. Eppendorf Centrifuge 5810R
14. Revco UltimaII -80°C Freezer
15. Thermolyne Locator 4 Liquid Nitrogen Freezer
16. Hemocytometer
17. Micropipet w/ P20 tips
18. Microscope

Growth Medium for Namalwa

RPMI 1640 medium with 2mM L-glutamine adjusted to contain 1.5g/L sodium bicarbonate, 4.5g/L glucose, 10mM HEPES, and 1.0mM sodium pyruvate (92.5%)
Characterized FBS (7.5%)
Pen-Strep (1X)

Cell Culture Procedure

A. Receipt of Frozen Cells and Starting Cell Culture

- 1) Immediately place frozen cells in liquid nitrogen storage until ready to culture.
- 2) When ready to start cell culture, quickly thaw the vial in a 37°C water bath.

- 3) As soon as ice crystals disappear, swab outside surface of the vial with 70% ethanol, then dispense contents of the vial into a 75cm² tissue culture flask and dilute with 40mL complete culture medium.
- 4) To incubate the culture, place the flask in a 37°C, 5% CO₂ humidified incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 2 million cells/mL.
- 2) Transfer 10mL of Namalwa culture into a new 75cm² tissue culture flask, and dilute with 30 mL of complete culture medium.
- 3) Record each subculture event as a passage.

C. Maintenance and Generation of Seed Stocks

- 1) Namalwa cells can be passaged every 3-4 days. Use 40mL of growth medium per 75cm² tissue culture flask.
- 2) Following first or second passage after receipt of cells and with sufficient number of cells to continue maintenance and expansion, the major portion of the flasks should be sub-cultured and a small portion should be set aside as a seed stock. The cell pellet (500 x g for 5 minutes) for the seed stock should be resuspended in freezing medium.
- 3) Cells in freezing medium are dispensed into cryovials (2 million cells per 1 mL aliquot) and frozen at -80°C in a Nalgene Cryo 1°C freezing container overnight.
- 4) Cryovials are transferred the next day to liquid nitrogen freezer for long-term storage.

D. Harvest

- 1) Passage cells until the desired number of cells is reached.
- 2) According to protocol described above under “Maintenance and Generation of Seed Stocks,” collect the cell pellet as for a seed stock (500 x g for 5 minutes).
- 3) Wash cell pellet by resuspension in PBS.
- 4) Examine viability using Trypan blue staining (SOP TP-7).
- 5) Collect cells by 500 x g 5 minutes centrifugation.

E. Sendai Virus infection

- 1) Start Namalwa culture as described above.
- 2) Thaw the stock Sendai Virus from the frozen vial (aliquoted upon receiving) in a 37°C water bath.
- 3) Start infection by adding the virus directly to the Namalwa culture flask to a final concentration of 200 HAU/mL. It will be ideal to infect cells at a density of 1 million cells/mL.
- 4) Continue to culture cells in a 37°C, 5% CO₂ humidified incubator for desired time.
- 5) Harvest cells by spinning at 500 x g for 5 minutes in Corning conical centrifuge tubes (15mL or 50mL).
- 6) Wash cell pellets by resuspension in PBS, and collect cells by 500 x g 5 minutes centrifugation.
- 7) Perform Percoll gradient purification according to **SOP: Cell purification using Percoll step gradients (04/01/2013)**.
- 8) Freeze purified Namalwa cells according to **SOP: Cryopreservation of hematopoietic cells from human leukapheresis product (02/07/2011)**.