

## Growth Protocol for L1210

Mouse L1210 lymphocytic leukemia cells (ATCC CCL219) and human leukemic MOLT-4 cells (ATCC CRL 1582) were grown in CO<sub>2</sub>-independent L-15 medium (Gibco BRL, Grand Island, NY, USA; Invitrogen catalog number 11415-064) supplemented with 2 gm/l dextrose and 10% fetal bovine serum (Gibco BRL).

## References

Grover NB, Eward KL, Helmstetter CE  
Synchronous cultures from the baby machine. A model for animal cells  
Mathematical Biosciences 190 (2004) 87-96